



## Can the L3 Linearity Solution in the HSQ kit be used for UV wavelength accuracy – FAQ

**The short answer is yes.**

Many other qualification vendors use only caffeine for wavelength **accuracy**, at 205nm and 273nm. Holmium oxide of course covers 241-641nm, and is directly traceable to NIST standards. However, if you want, there is no reason you could not simply take the spectrum of caffeine as well.

The L3 Solution provides a nice, mid-range absorbance value, compatible with most UV-Vis detectors. You could of course use any other concentration that gives a suitable response. The reference solution would be mobile phase, but there is very little difference if wanted wanted to use water.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-07-2024