





Sarcosine Analyzed with LCMS - AppNote


Separation of Potential Urine Biomarker from Isobaric β -Alanine

This developed LCMS method can separate Sarcosine from Beta- Alanine in serum and urine samples without using labor intensive sample derivatization. Since Sarcosine is considered a potential biomarker for prostate cancer risk and aggressiveness, it is essential to resolve and accurately quantify this compound in the presence of isobaric (same m/z) Beta-Alanine.

The developed method is Sensitive, Specific, Quantitative, and Reproducible ($\%RSD = 0.1$). It can be used in large scale studies with numerous samples (high throughput of the method due to simple sample preparation).



1. Sarcosine



2. β-Alanine

Peaks:



1. Sarcosine
2. β -Alanine

Method Conditions

Column: Cogent Diamond Hydride™, 4 μ m, 100Å

Catalog No.: 70000-15P-2

Dimensions: 2.1 x 15mm

Mobile Phase:

A: 50% Isopropyl Alcohol / 50% DI Water / 0.1% Acetic Acid

B: 97% Acetonitrile / 3% DI Water / 0.1% Acetic Acid

Gradient:

| Time (minutes) | %B |
|----------------|----|
| 0 | 75 |
| 3 | 75 |
| 4 | 65 |
| 5 | 65 |
| 10 | 20 |
| 12 | 75 |



Temperature: 50°C

Post Time: 5 minutes

Flow rate: 0.6mL / minute

Detection: ESI - POS - Agilent 6210 MSD TOF Mass Spectrometer

Injection vol.: 1µL

Sample Preparation: 10mg / L each of Sarcosine and Beta-Alanine in 50:50 A:B

***Note:** When Reversed Phase Columns were evaluated for their ability to separate Sarcosine from Beta-Alanine, both compounds eluted at the solvent front and were not separated. To achieve separation, a very intensive sample preparation has to be employed (e.g. derivatization) when using RP methods.*



Attachment

No 129 Sarcosine Analyzed with LCMS pdf 0.2 Mb [Download File](#)