


## Phosphorylated Sugars Analyzed with LCMS – AppNote

### Nucleotide Sugars UDP, ADP, TDP and CDP

UDP Hexanolamine was used as the Internal Standard in the analysis of Nucleotide Sugars in this rapid analysis Method. The Mobile Phase was designed to maximize the detector response in LCMS for maximum Chromatographic Efficiency. The simple “Inverse Gradient” which produces an Aqueous Normal Phase (ANP) HPLC Method, was required for the results shown.

This Method can be used in measuring Metabolite Concentration.



#### Peaks:

1. ADP-Glucose, *RT* = 3.68 min, monitored *MRM* transitions were *m/z* 588 to *m/z* 346
2. Proprietary Sugar Nucleotide, *RT* = 6.03 min, monitored *MRM* transitions were *m/z* 563 to *m/z* 321
3. Proprietary Sugar Nucleotide, *RT* = 6.20 min, monitored *MRM* transitions were *m/z* 606 to *m/z* 385
4. CDP-Glucose, *RT* = 7.13 min, monitored *MRM* transitions were *m/z* 564 to *m/z* 322
5. UDP Hexanolamine (internal standard), *RT* = 8.40 min, monitored *MRM* transitions were *m/z* 502 to *m/z* 258

*MRM* = Multiple Reaction Monitoring in LCMS/MS

### Method Conditions

**Column:** Cogent Diamond Hydride™, 4µm, 100Å

**Catalog No.:** 70000-15P-2

**Dimensions:** 2.1 x 150mm

#### Mobile Phase:

A: DI Water / 0.1% Ammonium Formate *pH* 7.2

B: 90:10 Acetonitrile / DI Water / 0.1% Ammonium Formate *pH* 6

#### Gradient:

Time (minutes)	%B
0	95
10.0	75

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

12.0	75
12.1	95
15.0	95

**Flow rate:** 0.3mL / minute

**Sample Preparation:** 400mL DI Water / 400mL Acetonitrile / 20mL of Stock Solution of each compound / 5mL of 12% ammonia

**Detection:** ESI - NEG - Agilent 6210 MSD TOF Mass Spectrometer

**Note:** Please note the addition of small amount of Ammonia to the Sample Matrix. The alkaline environment of the Sample Matrix assured Efficient and Symmetrical Peaks for all analytes.



## Attachment

**No 70 Phosphorylated Sugars Analyzed with LCMS pdf** 0.3 Mb [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)