


Carotenoids Analyzed with SFC – AppNote

Separation in Super Critical Fluid using a Bidentate C18 Column

Cogent TYPE-C Silica™ Stationary Phases with the unique hydrophobic, Silica Hydride Surface and a very wide range of Selectivities, can be used for Supercritical Fluid Chromatography (SFC).

In HPLC, Cogent TYPE-C Columns are known to retain many polar compounds in the Aqueous Normal Phase (ANP) Mode (*same mobile phases as HILIC*); in this Application Note the Separation behavior of the Bidentate C18 Column for three Carotenoids (Zeaxanthin, all trans Beta-carotene and 13 cis Beta-carotene) using SFC is shown to be similar.

High Efficiency and excellent Peak Symmetry were obtained with this Method when compared with others. An interesting result is the relatively long retention time of Zeaxanthin, which is eluted near all trans Beta-carotene. This unusual behavior displays the ability of Bidentate C18 to retain Polar Compounds more strongly than typical C18 phases in SFC.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation


9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Analytes



Peaks:

1. Zeaxanthin
2. All trans Beta-carotene
3. 13 cis Beta-carotene

Method Conditions

Column: Cogent Bidentate C18™, 4µm, 100Å

Catalog No.: 40018-25P

Dimensions: 4.6 x 250mm

Mobile Phase: 15% Methanol / 85% Carbon Dioxide

Outlet Pressure: 15MPa

Injection vol.: 5µL

Flow rate: 3mL / minute

Detection: UV @ 440nm

Notes: The small peaks in the chromatogram are impurities and / or degradation products in the samples. The Mobile Phase modifier is Methanol.



Attachment

No 43 Carotenoids Analyzed with SFC pdf 0.2 Mb [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com