



## Quinine Sulfate and Impurity Analysis with HPLC – AppNote

### API Separation from Impurity without Ion-Pair Reagents

Click [HERE](#) for Column Ordering Information.

The USP Method for Quinine Sulfate requires a Resolution of not less than 1.2 from its main impurity, Dihydroquinine. In the USP Method, the ion pair agents Methanesulfonic Acid and Diethylamine are used in the Mobile Phase. Ion pair agents are often needed to reduce peak tailing of basic analytes such as Quinine when conventional HPLC Columns are used.

In this Method, Quinine is separated from its main impurity with a Resolution of 2.6 using only 0.1% Trifluoroacetic Acid in the Mobile Phase.



#### Peaks:

1. Minor impurity
2. Quinine (API)
3. Dihydroquinine (Main impurity)

### Method Conditions

**Column:** Cogent Phenyl Hydride™, 4µm, 100Å

**Catalog No.:** 69020-7.5P

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Dimensions:** 4.6 x 75mm

**Mobile Phase:**

A: DI Water / 0.1% Trifluoroacetic Acid (TFA)

B: Acetonitrile / 0.1% Trifluoroacetic Acid (TFA)

**Gradient:**

Time (minutes)	%B
0	10
6	30
7	10

**Temperature:** 40°C

**Post Time:** 1 minute

**Injection vol.:** 10µL

**Flow rate:** 1.0mL / minute

**Detection:** UV @ 235nm

**Sample Preparation:**

*Stock Solution:* 1.0mg Quinine (90% by label claim) was dissolved in 1mL of 50:50 Solvent A / Solvent B mixture.

*Working Solution:* A 100µL aliquot of the stock was diluted with 900µL of 50:50 Solvent A / Solvent B mixture.

**t<sub>0</sub>:** 0.9 minutes

**Note:** Quinine is an alkaloid found naturally in the bark of Cinchona trees. It has been known for hundreds of years as an Antimalarial agent and a muscle relaxant. The Quechua inhabitants of Peru would grind the bark of Cinchona trees and mix it with sweetened water to reduce shivering in cold temperatures. Jesuits in Peru brought the bark to Europe where it was widely used to treat Malaria and became known as "Jesuit's bark."



**Attachment**

**No 154 Quinine Sulfate and Impurity Analysis with HPLC pdf** 0.3 Mb [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)