


## Water Soluble Vitamin Tablet Analysis with HPLC – AppNote

### Separation of Water Soluble Vitamins without Ion Pair Agents

This AppNote illustrates the Separation of four Hydrophilic, B Vitamins which were extracted from a commercially available tablet formulation. The data shows that these Water soluble B-Vitamins are well Separated from each other and from any Matrix Peaks. The main advantage of this Method over traditional Reversed Phase analyses for is that Ion Pair Agents are not needed to achieve acceptable results.

This Method can be used for LCMS applications. Furthermore, the Method's Precision is very good, as the overlay of 5 runs in the Figure below shows. Peak identities were confirmed by individual standards.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply


**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451


tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)


Website: [www.mtc-usa.com](http://www.mtc-usa.com)




1



2



3



4

### Peaks:

1. Niacinamide
2. Pyridoxine
3. Folic Acid
4. Cyanocobalamin

## Method Conditions

**Column:** Cogent Diamond Hydride™, 4µm, 100Å

**Catalog No.:** 70000-7.5P

**Dimensions:** 4.6 x 75mm

### Mobile Phase:

A: DI Water with 10mM Ammonium Formate

B: 95:5 Acetonitrile / Solvent A (v/v)

### Gradient:

Time (minutes)	%B
0	100
2	100
9	50
10	100

**Post time:** 3 minutes

**Injection vol.:** 1µL

**Flow rate:** 1.0mL / minute

**Detection:** UV @ 266nm

**Sample Preparation:** The Vitamin Tablet was ground and dissolved in 25mL of 50:50 10mM Ammonium Formate / Acetonitrile with 0.1% v/v 1N NaOH diluent. Solution was sonicated for 10 minutes and filtered through 0.45µm Nylon Syringe Filter (MicroSolv Tech Corp.).

**t<sub>0</sub>:** 1.0 minute

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

**Note:** The word “vitamin” was originally spelled “vitamine” when it was first coined by biochemist Casimir Funk. It was derived from the words “vital” and “amine” because it was believed at the time that all vitamins were chemical amines. The “e” was dropped from the word when it was discovered that this is not the case.



## Attachment

**No 171 Water Soluble Vitamin Tablet Analysis with HPLC pdf** 0.6 Mb [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)