

Before the first use of the Diamond Hydride column:

1. Condition it with 50:50 methanol / DI water for 30 minutes.

Conditioning steps - depends on your samples:

1. When using urine, plasma, or other biological samples etc:

A. Conditioning the column with 80% acetonitrile / 20% DI water/ 0.2% acetic acid (*or 0.1% formic acid*) is recommended for 20 minutes.

B. Cleaning with 50% MeOH / 50% DI water is recommended for every 10 injections for 5 minutes.

2. Using Standard Samples:

A: Conditioning with 80% acetonitrile / 20% DI water / 0.2% acetic acid (*or 0.1% formic acid*) is recommended for 20 minutes.

B. Cleaning with 50% MeOH / 50% DI water is recommended for every 100 injections or sooner if change in performance is noticed, for 5 minutes.

Storing the column:

1. For long term storage, 100% acetonitrile is best. Rinse with 10 column volumes.
2. For moderate length of time, 90:10 acetonitrile / DI water with 0.1% formic or acetic acid. Rinse with six column volumes.

Washing the column to remove adsorbed compounds:

The best method for washing a Diamond Hydride column is to use 50:50 Methanol / water with 0.1% formic acid as rinse solvent. Wash for 60 minutes using a slower flow rate than for your regular runs. Return to acetonitrile / DI water.

If the performance can still be improved, wash your column overnight using 50:50 isopropyl alcohol / DI water with 0.1% formic acid at a low flow rate.

If needed, uninstall the column and reinstall in a reverse direction and wash it with 50:50 methanol / DI water with 0.1% formic acid for 60 minutes.

Never let the column dry out.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com