

## Is the Cogent UDC Cholesterol HPLC column stable as a stationary phase – FAQ

**While cholesterol is a large molecule, the UDC-Cholesterol phase has proven to be stable** in 0.1% v/v TFA and 0.1% v/v H<sub>3</sub>PO<sub>4</sub> at acid pH's and up through pH 8.00 with 20mM NH<sub>4</sub>COOH **buffer** with NH<sub>4</sub>OH.

A pH in excess of pH 8.00 should be avoided as this phase has a mono dentate (*single point*) attachment to the silica matrix and an intermediate ester link in its structure.

*For higher pH's the Cogent Bidentate C18™ is recommended.*



Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-07-2024