


Superior Resolution, Reproducibility, & Peak Shapes Compared to USP Method

Click [HERE](#) for Column Ordering Information.

The USP assay method for Levothyroxine requires that a resolution of not less than 5.0 must be demonstrated between Levothyroxine and related compound Liothyronine. A chromatogram obtained from following the USP method using a Type-B Silica based L10 Column is shown in *Figure B*.

The average resolution between the two compounds over five runs is 2.8, which does not satisfy the system suitability for resolution for this assay. *Figure A* shows the five-run overlay obtained from a method developed with the Cogent Phenyl Hydride Column. The average resolution in this case was 5.3.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply


MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

**Peaks:**

1. Liothyronine Sodium
2. Levothyroxine Sodium

Method Conditions**Columns:**

Fig. A: Cogent Phenyl Hydride™, 4µm, 100Å

Fig. B: Type B Silica Based Column, 5µm, 100Å

Catalog Nos.:

Fig. A: 69020-7.5P

Fig. B: N/A

Dimensions:

Fig. A: 4.6 x 75 mm

Fig. B: 4.6 x 250 mm

Mobile Phase:

Fig. A:

A: DI Water / 0.1% Formic Acid (v/v)

B: 97% Acetonitrile / 3% DI Water / 0.1% Formic Acid (v/v)

Fig. B: 60% DI Water / 40% Acetonitrile / 0.05% Phosphoric Acid

Gradient:

Time (minutes)	%B
0	20
6	50
7	20

Temperature:

Fig. A: 35°C

Fig. B: Ambient

Injection vol.:

Fig. A: 2µL

Fig. B: 100µL

Flow rate:

Fig. A: 1.0 mL / minutes

Fig. B: 1.5 mL / minutes

Sample Preparation: Mix of Levothyroxine and Liothyronine standards.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Stock Solution: 0.4 mg Levothyroxine or Liothyronine dissolved with 1 mL 10 mM NaOH in 50:50 DI Water: Methanol.

Working Solution: Fig. A: Aliquots of stock solutions were mixed and diluted with 50:50 A:B to obtain concentrations of 40 mg / L and 4 mg / L for Levothyroxine and Liothyronine respectively.

Working Solution: Fig. B: Aliquots of stock solutions were mixed and diluted with the Mobile Phase to obtain concentrations of 10 mg /L and 0.2 mg / L for Levothyroxine and Liothyronine respectively.

Note: Levothyroxine is the L-isomer of the main thyroid hormone Thyroxine (T4). It is used as a replacement for the Thyroxine that is deficient in patients with hypothyroidism. Liothyronine is the L-isomer of another thyroid hormone, Triiodothyronine (T3). T3 is produced from T4 and is the metabolically active form of the hormone.



Attachment

No 128 Levothyroxine Assay pdf 0.3 Mb [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com