



A: DI Water / 10mM Sodium Phosphate Monobasic

B: 80% Acetonitrile / 20% DI Water / 10mM Sodium Phosphate Monobasic

**Gradient:**

Time (minutes)	%B
0	0
15	100

**Temperature:** 40°C

**Injection vol.:** 10µL

**Flow rate:** 1.5mL / minute

**Detection:** UV @ 280nm

**Notes:** Benzonatate is a non-narcotic oral antitussive (cough suppressant) drug, taken orally, which works by anesthetizing the tissues of the lungs and pleura responsible for the cough reflex. It is sold under the brand name Tessalon Perles (gelcaps) or Tessalon Capsules.



**Attachment**

**No 80 Benzonatate Excipients & Degradants pdf** 0.1 Mb [Download File](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)