



Alprazolam Analyzed with HPLC - AppNote

Robust Separation of API from USP Internal Standard

The USP Assay Method for Alprazolam uses a bare Silica Column and a complex Mobile Phase consisting of Acetonitrile, Chloroform, Butyl Alcohol, and Acetic Acid. In this Method, a simple LCMS compatible Mobile Phase is used and produces excellent Peak shapes for both the API and its USP internal standard.

Furthermore, a resolution of 4.3 was obtained between the two Peaks, which meets the USP system suitability of $R_s \geq 2.0$. Two impurity Peaks are also observed, which further illustrates the resolution capabilities of the analysis. This Method illustrates how the MS-compatible HPLC-UV Methods described in various application notes can be successfully adapted for LCMS.





Peak 1

Peak 2

Peaks:

1. Triazolam (internal standard)
2. Alprazolam (API)
- 3, 4. Impurities

Method Conditions

Column: Cogent Diamond Hydride™, 4μm, 100Å

Catalog No.: 70000-7.5P

Dimensions: 4.6 x 75mm

Mobile Phase:

A: DI Water / 0.1% Formic Acid (v/v)

B: Acetonitrile / 0.1% Formic Acid (v/v)

Gradient:

Time (minutes)

%B



0	95
1	95
6	50
7	95

Post Time: 3 minutes

Injection vol.: 1 μ L

Flow rate: 1.0mL / minute

Detection: UV @ 254nm

Sample Preparation:

Tablet: A 0.25mg strength generic Xanax® tablet was ground and added to a 10mL volumetric flask. After diluting with 50% Solvent A / 50% Solvent B, it was sonicated for 10 minutes. A portion was filtered with a 0.45 μ m Nylon Syringe Filter (MicroSolv Tech Corp.).

Internal Standard: 1mg / mL Triazolam in Methanol diluent.

Working Solutions: 20 μ L of the internal standard and 980 μ L of the tablet extract were mixed. Peak identities were confirmed by individual solutions of the tablet extract and the internal standard.

t_0 : 0.9 minutes

Note: Alprazolam is a member of the benzo-diazepine class of compounds, prescribed to treat a variety of anxiety-related conditions.



Attachment

No 183 Alprazolam Xanax Analyzed with HPLC pdf 0.6 Mb [Download File](#)