



## What are the dimensions and specifications of the 10L reservoir bottles - FAQ

### **Our 10L (10,000ml) safety coated solvent reservoir or mobile phase storage bottles have the following specifications:**

*Material: Borosilicate Glass Code 7740*

*Safety Coating: Proprietary epoxy plastic*

**Capacity:** *10,000 mL*

*Bottle Style: Storage*

*Neck Style: Wide Mouth*

*Top Style: Threaded*

*Cap Style: Plug Seal*

*Cap Material: Polypropylene*

*Bottle Shape: Cylindrical*

*Cap Color: Orange*

*Thread Finish: 45 GL*

*Graduation Range: 2000-9000 mL (approx)*

*Graduation Interval: 1000 mL (approx)*

*Bottle Feature: Graduated*

*Feature: Reusable*

*Height: 410 mm (approx)*

*Outer Diameter: 234 mm (approx)*

*Note: Bottles larger than 2L should NOT be used with bottle top filter units, or in other applications involving vacuum pressure, as breakage may occur. DO NOT tighten caps immediately after autoclaving as the vacuum resulting from cooling can cause breakage. A vented cap is NOT recommended for use on the 10L bottle for venting during autoclaving.*

[HPLC Reservoir Bottle Ordering Information](#)

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 09-19-2024