



## If I use a buffer at pH 9 at 30 to 35C with the Zero Flow capillaries what lifetime should I expect from them - FAQ

**If you use a pH buffer even with moderate ionic strength, at pH 9, the capillaries would *not* be expected to last very long.**

It is hard to predict how many runs you get because it will depend on your run times, your conditioning step, etc. However, these capillaries are not designed for a pH higher than 8. The silica was and the organo-silane bond of the polyacrylamide is susceptible to hydrolysis, degradation and pitting at high ionic strength buffers above pH 8.

Click [HERE](#) for capillaries for CE ordering information and pictures

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)