

## Calculate the internal volume of tubing – How to

To calculate the internal volume (V), you need to know the length of tubing (L) and the inner diameter (ID). You then use the formula for volume of a cylinder.

Convert L and ID into cm first. (1 inch = 2.54 cm). This gives V in cm<sup>3</sup>. 1 cm<sup>3</sup> = 1 mL. You can then convert into µL if desired by 1000µL = 1mL.

$$V = \pi (ID/2)^2 L$$



[Tubing Product Page](#)

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-05-2024