



## Partisil C18 stationary phase names: ODS1 ODS2 or ODS3

Sometimes, you may see a Partisil HPLC column referred to as “ODS.” This means the column has an Octadecyl silane (C18) stationary phase. In addition, you may see ODS followed by a number designation:

**Partisil™ ODS** is a low carbon load C18 phase, with intentionally unbonded surface silanols adding to the selectivity.

**Partisil™ ODS2** is a relatively high carbon load C18 bonded silica phase which makes it the most non-polar and, therefore, the most retentive of the Partisil reversed-phase columns. The high sample load **capacity** is suitable for preparative work.

**Partisil™ ODS3** is an intermediate carbon load C18 phase used for pharmaceutical, natural products, food, biological and environmental pollutants applications.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 04-05-2024