



Topiramate, Retention of a Sulfonate Drug by LCMS – AppNote

Retention of the Sulfonate Drug Topiramate

The structure of Topiramate leads us to think due to the Amino group and negative Log P value of -0.8, it will not be retained in Reversed Phase HPLC. However, in this Method, the compound is well retained.



Peak:

Topiramate 340.1062 m/z (M + H)⁺

Method Conditions:

Column: Cogent Bidentate C18™, 4μm, 100Å

Catalog No.: 40018 -05P-2

Dimensions: 2.1 x 50mm

Mobile Phase:

A: DI Water / 0.1% Formic Acid

B: Acetonitrile / 0.1% Formic Acid

Gradient:

Time (minutes)	%B
0	10
1	10
4	90
5	90
6	10
7	10

Injection vol.: 1µL

Flow rate: 0.3mL / minute

Detection: ESI - pos - Agilent 6210 MSD TOF Mass Spectrometer

Sample Preparation: Topiramate 0.01mg / mL in 50% A / 50% B solvent mixture.

Notes: *Migraines are a common, disabling neurological disorder and are often accompanied by one or more of the following disabling symptoms: visual disturbances, nausea, vomiting, dizziness, extreme sensitivity to sound, light, touch and smell, and tingling or numbness in the extremities or face. Migraines involve a complex interchange between various brain regions, including the hypothalamus and brainstem nuclei that modulate pain signaling. The headache phase involves activation of the Trigemino-vascular system. Topiramate modulates the Trigemino-vascular signaling that is effective in migraine prevention. It is also used as an Antiepileptic drug (AED) and is structurally distinct from other AEDs as it is derived from D-fructose, a naturally occurring sugar moiety, and has sulfamate functionality*



Attachment

A377 Retention of a Sulfonate Drug by LCMS pdf 0.1 Mb [Download File](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 03-05-2024