

Why use Ammonium Acetate in ANP and LCMS – Tips & Suggestions

Ammonium acetate can be a very useful additive for HPLC methods.

With ammonium acetate, a near-neutral pH is obtained with this additive, so it is a popular choice for keeping compounds such as most organic acids ionized, which is generally preferred for Aqueous Normal Phase (ANP) HPLC methods.

In addition, it is LC-MS compatible, in contrast to a near neutral pH phosphate buffer. However, it also is thought to adsorb to silica based stationary phase materials and therefore the column should thereafter be dedicated to methods using ammonium acetate.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 09-26-2024