

## What is 0.1 % w/v ammonium acetate in terms of molarity - FAQ

Percentage given as w/v means **weight over volume**. This is typically used for solids dissolved in a liquid.

Ammonium acetate is a solid and is often used in HPLC mobile phases. Suppose you have a 0.1 % (w/v) ammonium acetate solution but your SOP requires you to state the concentration in mM. How do you convert? You will need the molecular weight of the solute (77.1 g/mol in this case.)

$$0.1 \% (w/v) = \frac{0.001 \text{ g}}{\text{mL}}$$

$$\left[ \frac{0.001 \text{ g}}{\text{mL}} \right] \left[ \frac{\text{mol}}{77.1 \text{ g}} \right] \left[ \frac{1000 \text{ mL}}{\text{L}} \right] = \frac{0.0130 \text{ mol}}{\text{L}} = 13 \text{ mM}$$

1. Convert the percent to grams per mL
2. Change grams to mol (use the molecular weight)
3. Change mL to L
4. You will get mol per liter, or M. This may be expressed as mM by multiplying by 1000.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 07-22-2024