

Suggestions on how to separate capsaicinoids in hot sauce by HPLC - Tips & Suggestions

If you have a hot sauce extract containing nordihydrocapsaicin, capsaicin, dihydrocapsaicin, and homodihydrocapsaicin, the Cogent Bidentate C18 2.0™ HPLC column is a good choice for sufficient separation when the compounds are analyzed by LC-MS; look for the $[M+H]^+$ ions.

The following gradient should be a good starting point:


Mobile Phase: A: DI H₂O / 0.1% formic acid (v/v)

B: Methanol

Gradient:

time (min.)	%B
0	10
1	10
9	70
10	70
11	10

You may get EICs which look like this:



Adapted from the Chrom Resource Center

Copyright 2024, All Rights Apply

1. nordihydrocapsaicin, 294.2064 m/z $[M+H]^+$ **MicroSolv Technology Corporation**

2. capsaicin, 306.2064 m/z $[M+H]^+$ 9158 Industrial Blvd. NE, Leland, NC 28451

3. dihydrocapsaicin, 308.2220 m/z $[M+H]^+$ tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

4. homodihydrocapsaicin, 322.2377 m/z [M+H]⁺

For peaks 1 and 2, you will be able to obtain separate peaks for each compound in the EICs so there will not be an issue of interference from co-elution.

In terms of sample preparation, try using a Soxhlet extraction procedure.

Click [HERE](#) for Cogent Bidentate C18 HPLC Column Ordering Information.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com