



## Is formic acid compatible with borosilicate glass – FAQ

**Yes, formic acid is generally compatible with borosilicate glass.**

There are two aspects that one needs to be mindful of when storing formic acid in glass vials or glass inserts:

1. Formic acid gradually decomposes to carbon monoxide (CO) and water, building up pressure inside the container. The warmer the solution the faster the process. Therefore, formic acid is typically supplied in plastic bottles and is usually refrigerated.
2. Universal to all solutions stored in glass, is leaching alkali metals, mainly sodium, boron and silicic acids from the surface of the glass containment into the stored solution.

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)