



A reproducible assay for detection of a non-ionic surfactant.

Triton X-100 is used in the sample purification steps of various materials as well as for cleaning apparatus in pharmaceutical manufacturing. Detection and quantification are very important to identify any residual amounts of this compound.

The overlay of 5 injections (below) demonstrates reproducibility in separating the polyoxyethylene mixture.



Peak:

Triton X-100

Method Conditions:

Column: Cogent Bidentate C8™, 4µm, 100Å

Catalog No.: 40008-10P

Dimensions: 4.6 x 100mm

Mobile Phase:

A: DI water with 0.1% formic acid

B: Acetonitrile with 0.1% formic acid

Time (minutes)	%B
0	98
2	95
3	95
5	95
10	98

Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Injection vol.: 1µL

Flow rate: 1.0mL / minute

Detection: UV @ 230nm

Sample Preparation: Triton X-100, 10uL / mL in DI water

Note: Triton X-100 is one of the most widely used nonionic surfactants for lysing cells to extract protein and other cellular organelles. Triton X-100 ($C_{14}H_{22}O(C_2H_4O)_n$) has a hydrophilic polyethylene oxide chain and an aromatic hydrocarbon lipophilic or hydrophobic group. The hydrocarbon group is a 4-phenyl group.



Printed from the Chrom Resource Center

Copyright 2024, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com