




## Organic Acids An Important Class of Compounds







**Notes:** Chromatogram presented is adapted from: J.J.Pesek, M.T. Matyska, S. M. Fisher, T. R. Sana, Journal of Chromatography A, 1204 (2008) p53.

## **Method Conditions**

## Column: Cogent Diamond Hydride™, 4µm, 100Å

Catalog No.: 70000-15P-2

Dimensions: 2.1 x 150 mm

Solvents: A: DI H<sub>2</sub>O + 10 mM ammonium formate B: 90% acetonitrile/ 10% DI H<sub>2</sub>O/ 10 mM ammonium formate

dient:	time (min.)	%B
	0-3	90
	3-6	70
	6-7	70
	7-7.1	30
	7.1-8	90

Post Time: 5 min

Grad

Flow rate: 0.4 mL/min

Detection: ESI - neg - Agilent 6210 MSD TOF mass spectrometer

Sample: Sample stock solutions were made in acidified DI H<sub>2</sub>O at a concentration of 0.2 - 0.3 mg/mL. Samples for analysis were made by diluting the stock 1:100 in 50:50 solution A and B.

Peaks: 1. Maleic Acid, 115 m/z (M-H)<sup>-</sup>

- 2. Trans Aconitic Acid, 173 m/z (M-H)<sup>-</sup>
- 3. Aconitic Acid, 173 m/z (M-H)
- 4. Impurity 173 m/z (M-H)<sup>-</sup>
- 5. Fumaric Acid, 115 m/z ((M-H)<sup>-</sup>
- 6. Citric Acid, 191 m/z (M-H)-
- 7. Oxaloacetic Acid, 131 m/z (M-H)<sup>-</sup>

## Discussion

Organic acids, an important class of metabolites, are very well retained and separated in this method. Figure A shows the combined EICs for five small organic acids: maleic, aconitic, fumaric, citric and oxaloacetic acids. Two of these acids are isobaric compounds; maleic and fumaric acids (115 m/z) and are separated very well. There are three peaks at 173 m/z. The peaks were assigned to trans (peak 2) and cis (peak 3) isomers of the aconitic acid and an impurity (peak 4). Citric and oxaloacetic acids were not resolved using the developed mobile phase and gradient, but can be very easy assigned using MS detection. Figure B shows the EIC at 131 m/z for oxaloacetic acid for integration purpose.

APP-A-100



9158 Industrial Blvd NE Leland, NC 28451 p: 1.732.380.8900 f: 1.910.769.9435