



Nucynta® (Tapentadol HCI) Tablet

Assay method for analgesic amine





Note: Tapentadol is an analgesic compound used to treat moderate to severe pain. It efficacy is due to two modes of action: One is an agonist of the μ -opioid receptor and another as a norepinephrine reuptake inhibitor.



Column: Cogent Diamond Hydride™, 4µm, 100Å

Catalog No.: 70000-7.5P

Dimensions: 4.6 x 75 mm

Mobile Phase: A: DI H₂O / 0.1% formic acid (v/v) B: Acetonitrile / 0.1% formic acid (v/v)

Gradient:	time (min.)	%B
	0	95
	1	95
	6	40
	7	95

Post Time: 3 min

Injection vol.: 1µL

Flow rate: 1.0 mL/min

Detection: UV 271 nm

Sample: 75mg strength Nucynta* tablet was ground and weighed in a 25mL volumetric flask. A portion of 50/50 solvent A/solvent B diluent was added and the flask was sonicated 10 min. It was then diluted to mark and filtered with a 0.45µm nylon syringe filter (MicroSolv Tech Corp.). The filtrate was diluted 1:5 for HPLC injections.

Peak: 1. Tapentadol HCI

t₀: 0.9 min

Discussion

Tapentadol can be a problematic compound for HPLC analysis due to the amine. Tertiary amines are often particularly difficult to obtain a good peak shape for. Using reversed phase methods, notable tailing is observed in several published reports in the literature. The Cogent Diamond Hydride column, however, obtains a sharp peak for this analyte due to its unique retention mode.

Data from two column lots shown in the figure illustrates the reproducibility of the material from batch to batch.

MANUFACTURED BY: MICROS UV TECHNOLOGY CORPORATION

9158 Industrial Blvd NE Leland, NC 28451