



## Effexor® (Venlafaxine) Capsule

## LC-MS compatible method for phenethylamine compound





**Note:** Venlafaxine is a serotonin-norepinephrine reuptake inhibitor used to treat various depressive and anxiety disorders. It is currently marketed by Pfizer as Effexor®. It is a phenethylamine and shares structural similarities with other compounds in this class, such as amphetamine, methamphetamine, and MDMA.

## **Method Conditions**

Column: Cogent Diamond Hydride™, 4µm, 100Å

**Catalog No.:** 70000-7.5P **Dimensions:** 4.6 x 75 mm

**Mobile Phase:** A: DI  $H_2O$  / 0.1% formic acid (v/v)

B: Acetonitrile / 0.1% formic acid (v/v)

 Gradient:
 time (min.)
 %B

 0
 95

 1
 95

 6
 50

 7
 95

Post Time: 3 min
Injection vol.: 1µL
Flow rate: 1.0 mL/min
Detection: UV 226 nm

Sample: 75mg strength Effexor® Extended Release capsule contents were added to a 25mL volumetric flask. A portion of 50/50 solvent A/solvent B diluent was added and the flask was sonicated 10 min. It was then diluted to mark and mixed. A portion was filtered with a 0.45µm nylon syringe filter (MicroSolv Tech Corp.) and diluted 1:50.

Peak: 1. Venlafaxine

to: 0.9 min

## **Discussion**

The USP assay method for venlafaxine capsules uses triethylamine and phosphoric acid in the mobile phase, both of which are incompatible with LC-MS. The system suitability for venlafaxine tailing factor is 2.0, indicating the compound has a tendency for tailing. Here a sharp symmetrical peak is observed using formic acid. Data from two column lots is shown to illustrate reproducibility.