



Tizanidine HCl Method Transfer

Use of near UHPLC column for improved results



Tizanidine

Note: Tizanidine is a centrally acting α_2 -adren-
ergic agonist used to treat spasms, cramping,
tightness of muscles, and related conditions. It is
available under the trade name Zanaflex® as well
as generic versions.

Method Conditions

Column: Fig. A: Cogent Diamond Hydride 2.0™, 2.2µm, 120Å

Fig. B: Cogent Diamond Hydride™, 4µm, 100Å

Catalog No.: Fig. A: 70200-05P-2; Fig. B: 70000-05P-2

Dimensions: 2.1 x 50 mm

Mobile Phase: A: DI H₂O / 0.1% formic acid (v/v)

B: Acetonitrile / 0.1% formic acid (v/v)

Gradient:	time (min.)	%B
	0	95
	1	95
	6	40
	7	95

Post time: 3 min

Injection vol.: 0.2 microL

Flow rate: 0.3mL/min

Detection: UV 230 nm

Sample: 4mg strength tizanidine HCl tablet was ground and weighed
in a 10mL volumetric flask. A portion of 50/50 solvent A/
solvent B diluent was added and the flask was sonicated 10
min. It was then diluted to mark and filtered with a 0.45µm
nylon syringe filter (MicroSolv Tech Corp.).

Peak: 1. Tizanidine HCl

t₀: 0.6 min

Discussion

This challenging compound has several amine groups and can be
problematic to analyze by HPLC. Here an excellent peak shape
is obtained using the Cogent Diamond Hydride 2.0 column. The
efficiency is higher on the 2.2µm phase compared to a standard 4µm
column, leading to increased sensitivity. The method conditions are
compatible with LC-MS as well.