


ANP

Sugars, Sugar Alcohols & Amino Acids

Glucose, Mannitol, Aspartic Acid, Alanine and Arginine



| 102 | 1. | .2 | 3.4 | F |
|------|----|----|-----|---|
| 0.9- | | - | | 5 |
| 0.8- | | | | |
| 0.7- | | | | |
| 0.6- | | | | |
| 0.5- | | | | |
| 0.4- | | | | |
| 0.3- | | | | |
| 0.2- | | | | |
| 0.1- | | | | |

Notes: Sugar alcohols differ from sugars in that the carbonyl group of the former has been reduced to an alcohol group. They don't raise blood sugar levels to the same extent that regular sugars do and hence are often used in diabetic foods when a sweet taste is still desired. Amino acids have many biological functions and are the building blocks of peptides and proteins.

Method Conditions

Column: Cogent Diamond Hydride™, 4µm, 100Å

Catalog No.: 70000-15P-2

Dimensions: 2.1 x 150 mm

Solvents: A: 50% DI $H_2O/$ 50% acetonitrile/ 0.1% formic acid/ 100µM sodium acetate B: Acetonitrile/ 0.1% formic acid

| Gradient: | time (min.) | %B |
|-----------|-------------|----|
| | 0 | 90 |
| | 1 | 90 |
| | 3 | 80 |
| | 4 | 80 |
| | 5 | 50 |
| | 6 | 50 |
| | 7 | 25 |
| | 10 | 25 |
| | 11 | 90 |

Post Time: 5 min

Injection vol.: 1µL

Flow rate: 0.6 mL/min

Detection: ESI - POS - Agilent 6210 MSD TOF mass spectrometer

Peaks: 1. Glucose

- 2. Mannitol
- 3. Aspartic Acid
- 4. Alanine
- 5. Arginine

Discussion

Metabolomics analysis encompasses a vast array of molecules with disparate hydrophobic/hydrophilic properties. Many of these analytes are in fact quite polar and problematic to retain by conventional reversed phase chromatography. However, the Cogent Diamond Hydride column is quite well-suited to their separation. This application demonstrates this property through the separation of five notable polar compounds that may be of interest to metabolomics separations.

APP-A-366

MICROS IV TECHNOLOGY CORPORATION

9158 Industrial Blvd NE Leland, NC 28451 p: 1.732.380.8900 f: 1.910.769.9435 customers@mtc-usa.com www.Cogent-HPLC.com