

Calculating theoretical plates N at half height - How To

The theoretical plate number $N_{0.5}$ can be calculated from a chromatogram using the equation below:

$$N_{0.5} = 5.54 (t_R/w_{0.5})^2$$

Where t_R is retention time and $w_{0.5}$ is the peak width at half its height. Remember to use the same units (minutes or seconds) for both values (t_R and $w_{0.5}$).



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com