



Why use Ammonium Acetate in ANP and LCMS - Tips & Suggestions

Ammonium acetate can be a very useful additive for HPLC methods.

With ammonium acetate, a near-neutral pH is obtained with this additive, so it is a popular choice for keeping compounds such as most organic acids ionized, which is generally preferred for Aqueous Normal Phase (ANP) HPLC methods.

In addition, it is LC-MS compatible, in contrast to a near neutral pH phosphate buffer. However, it also is thought to adsorb to silica based stationary phase materials and therefore the column should thereafter be dedicated to methods using ammonium acetate.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com