



## Coricidin HBP Night® Tablet - AppNote

### Acetaminophen, Dextromethorphan, and Chlorpheniramine

Click [HERE](#) for Column Ordering Information.

This Method demonstrates the potential of the Cogent Phenyl Hydride 2.o™ Column for analysis of amine-containing analytes in a real formulation. Chlorpheniramine and Dextromethorphan can often be problematic in terms of peak tailing due to the tertiary amines, but this column and method produces excellent peak shapes for both compounds in this tablet formulation.

A high UV wavelength of 310 nm was chosen at the beginning of the gradient in order to reduce the Acetaminophen peak height, since this analyte is present in much higher concentration than the other two.



#### Peaks:

1. Acetaminophen
2. Chlorpheniramine
3. Dextromethorphan

### Method Conditions

**Column:** Cogent Phenyl Hydride 2.o™, 2.2μm, 120Å

**Catalog No.:** [69220-05P-2](#)

**Dimensions:** 2.1 x 50 mm

**Mobile Phase:**

A: DI Water / 0.1% Trifluoroacetic Acid (TFA) (v/v)

B: Acetonitrile / 0.1% Trifluoroacetic Acid (TFA) (v/v)

Time ( <i>minutes</i> )	%B
0	5
2	5
9	80
11	80
12	5

**Post Time:** 5 minutes

**Injection vol.:** 1µL

**Flow rate:** 0.3mL / minute

**Detection:** UV @ 310nm (0-2 min), @ 265 nm (2-12 min)

**t<sub>0</sub>:** 0.7 minutes

**Samples:** Coricidin High Blood Pressure (HBP) Night® tablet containing 500mg Acetaminophen, 15mg Dextromethorphan HBr, and 2mg Chlorpheniramine maleate was ground and added to a 50mL volumetric flask. A portion of 50/50 solvent A/solvent B was added and the flask was sonicated for 10 minutes. Solution was then diluted to mark and mixed well. A portion was filtered with a 0.45µm Nylon Syringe Filter (MicroSolv Tech Corp.). This filtrate was used for injections and peak identities were confirmed with individual standards.



**Attachment**

**No 324 Coricidin HBP Night Tablet pdf** 0.2 Mb [Download File](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)