



What is the structure of the Cogent Phenyl Hydride bonded ligand - FAQ

The Cogent Phenyl Hydride phase is a 4 carbon chain with a phenyl group at the end attached directly to the silica hydride surface with a single, direct silicon-carbon bond.

With bonds in traditional silica B based materials, the organo-silane starting material is Si-O-Si-C4-Phenyl and this reagent is reacted with the silica to form the final stationary phase which includes a Si-O-Si-C4-Phenyl ligand.

This is in contrast with the TYPE-C phases in that the Si-C bond is formed directly as part of the reaction resulting in Si-C4-Phenyl without the oxygen in the link which is very susceptible to hydrolysis.

[More information on Phenyl Hydride HPLC Columns](#)



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com