




# Glyburide and Metformin

Separation of “Highly Polar” and “Non-Polar” Compounds in one isocratic run



**Notes:** Elution order was confirmed by LC-MS, APCI+, with single ion monitoring Metformin (m/z 130) and Glyburide (m/z 369).

## Method Conditions

**Column:** Cogent Bidentate C18™, 4µm, 100Å

**Catalog No.:** 40018-75P

**Dimensions:** 4.6 x 75 mm

**Mobile Phase:** A: 15% DI H<sub>2</sub>O/ 85% acetonitrile/ 0.5% formic acid

B: 50% DI H<sub>2</sub>O/ 50% acetonitrile/ 0.5% formic acid

**Injection vol.:** 1µL

**Flow rate:** 0.5 mL/min

**Detection:** UV 254 nm

**Sample: Stock Solution:** 100µg/µL glyburide and metformin

## Discussion

The polar compound, Metformin, and the nonpolar compound Glyburide, can be retained on the same stationary phase (see A & B), utilizing the unique properties of the hydride based Cogent Bidentate C18 column. Depending on the mobile phase composition either Metformin or Glyburide can be retained longer.