



Zinc-EDTA Complex

No ion pair agents needed



Note: EDTA (ethylenediaminetetracetic acid) has been used to promote metal uptake from contaminated soil. As such, the ESI-MS technique described here can be used for the determination of Zn-EDTA complexes in soil samples.

Method Conditions

Column: Cogent Diamond Hydride™, 4µm, 100Å

Catalog No.: 70000-15P-2

Dimensions: 2.1 x 150 mm

Mobile Phase: A: DI H₂O / 0.1% formic acid (v/v)
B: Acetonitrile / 0.1% formic acid (v/v)

Gradient:	time (min.)	%B
	0	90
	5	20
	8	20
	9	90

Post time: 2 min

Injection vol.: 1µL

Flow rate: 0.40 mL/min

Detection: ESI - NEG - Agilent 6210 MSD TOF mass spectrometer

Sample: A soil sample was spiked with Zn-EDTA complex at a level of 2000µm. After extraction with DI H₂O (shaking for 24 hours), the sample was filtered using a 0.45µm syringe filter (MicroSolv Tech Corp.) and diluted with acetonitrile 1:10 before injection.

Peak: Zn-EDTA complex 354.7 m/z

t₀: 0.9 min

Discussion

Using conventional analytical methods, retention of metal-EDTA complexes is accomplished using ion pair reversed phase chromatography. However, the ion pair agents used in the mobile phase are not compatible with mass spectrometry. In this LC-MS method using the Cogent Diamond Hydride column, only formic acid is needed in the mobile phase in order to obtain retention of a Zinc-EDTA complex. The figure shows an EIC of the analyte spiked in a soil extract matrix.