

## Can BASIK brand screw caps be used in ion chromatography – FAQ

While BASIK™ brand screw caps are not designed for ion chromatography, they may work for some applications. But I would not use them sensitive methods.

We do suggest it might be worthwhile for maximum performance in IC to use our ultra-pure polyethylene screw caps, see image below.



*Note that these caps, are a single use cap and will not reseal after injection.*

Click [HERE](#) for more information about IC compatible caps.

Click [HERE](#) for PP vials that are IC compatible.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-08-2024