



## How similar fittings for chromatography can differ from different brands in HPLC – Tech Information

Of course, the quality of the manufacturing is Important but besides that the most important distinguishing features are:

1. Geometry of the receiving port you are using.
2. Threads call out on the fitting and receiving port, count per inch and total length of threads.
3. Material the fitting is made out of for purposes of chemical inertness, temperature the fitting will be used for.
4. Pressure the fitting will be required to hold.
5. Size of the tubing the fitting will be holding.

Click [HERE](#) for More Information on Fittings for chromatography

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-05-2024