



## What is the maximum temperature that the Endure brand of fittings can be used - FAQ

150°C for Endure™ carbon PEEK fittings and 100°C for most other PEEK tubing connector fittings also known as male nuts.

Click [HERE](#) for ordering information for tubing connectors

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-05-2024