



What is a Zero Flow CE capillary mostly used for – FAQ

There are many uses for Zero Flow™ capillary in CZE. One use is to allow analytes to separate based only on electrophoretic mobility. The Zero Flow™ capillary is an excellent choice for this technique because, regardless of the pH you use to ionize your analytes, the flow will still be zero.

Another excellent use is cIEF.

Click [HERE](#) for Zero Flow and Controlled Flow capillaries ordering information and pictures.

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-18-2024