



What is the advantage of using a Pharmafix type of septa for headspace analysis – FAQ

The advantage is that you get a better seal with the “Pharmafix” seal. However, you should be aware that it may not be as inert as a silicone/PTFE septa designed for headspace analysis.

Click [HERE](#) for ordering information and pictures of headspace caps

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-05-2024