



## Are there any solvents that I cannot use with the knitted reactor coils KRC and knitted delay tubes KDT - FAQ

Yes, you cannot use any perfluorinated solvents or any solvent not compatible with PTFE with the knitted reactor coils KRC and knitted delay tubes KDT.

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)