



For Normal Phase HPLC which syringe filter do you recommend – FAQ

PTFE is the best membrane type for use with non polar, organic solvents required for Normal Phase Chromatography NP. Click [HERE](#) for AQ™ brand PTFE syringe filter ordering information and pictures.

Attachment

MicroSolv filters equivalency study pdf [Download File](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 06-05-2024