



Using one syringe filter for many samples will I get sample carry over – FAQ

You will most likely get some kind of carry over. It is not recommended to use one filter for more than one sample.

[Click HERE for AQ syringe filter ordering information and pictures](#)

Attachment

MicroSolv filters equivalency study pdf [Download File](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-13-2024