



## What is the difference between the Beta Tools and other similar pliers on the market – FAQ

The Beta Tools™ are pairs of **precision** pliers not unlike many others. But what makes them very different from others who try to mimic it are the center holes. These **precision** drilled holes are very difficult to create each and every time as they have to be “dead center” for it to work properly.

Other manufacturers offer similar **precision** pliers with a drilled hole for tubing straightening but their tolerances of “dead center” are often not satisfactory. The Beta Tool™ has been in production for over 25 years and the manufacturing and QC process make it the tubing straightener of choice of labs and service technicians worldwide.

[Ordering information for the Beta Tool2™](#)

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)