



What are the dimensions and specifications of the 10L reservoir bottles – FAQ

Our 10L (10,000ml) safety coated solvent reservoir or mobile phase storage bottles have the following specifications:

Material: Borosilicate Glass Code 7740

Safety Coating: Proprietary epoxy plastic

Capacity: *10,000 mL*

Bottle Style: Storage

Neck Style: Wide Mouth

Top Style: Threaded

Cap Style: Plug Seal

Cap Material: Polypropylene

Bottle Shape: Cylindrical

Cap Color: Orange

Thread Finish: 45 GL

Graduation Range: 2000-9000 mL (approx)

Graduation Interval: 1000 mL (approx)

Bottle Feature: Graduated

Feature: Reusable

Height: 410 mm (approx)

Outer Diameter: 234 mm (approx)

Note: Bottles larger than 2L should NOT be used with bottle top filter units, or in other applications involving vacuum pressure, as breakage may occur. DO NOT tighten caps immediately after autoclaving as the vacuum resulting from cooling can cause breakage. A vented cap is NOT recommended for use on the 10L bottle for venting during autoclaving.

[*HPLC Reservoir Bottle Ordering Information*](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-18-2024