



## Can I autoclave the plastic Chrom Syringes -FAQ

**Yes, the Chrom Syringes are fully compatible in an autoclaving application.**

If you are using the chromatography syringes with syringe filters, please verify that the filter is also suitable for the autoclaving conditions as well.

[Luer-Lock Syringe Product Page](#)

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-15-2024