

What is the largest ID for PEEK columns that you offer – FAQ

The largest I.D. available is 7.5mm.

The corresponding lengths offered for a column of this ID are 30, 50, 100, 150, 250, 300 mm. All of these columns are 5/8" O.D.



[PEEK Hardware Product Page](#)

Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 06-05-2024