

## Can I use the tubing plier cutters with 1/16 or 1/8 OD stainless steel HPLC tubing - FAQ

This tool is compatible with both 1/16 inch and 1/8 inch OD stainless steel tubing. The tool can be used to score stainless steel tubing and then you can manually break the tubing at the score. This tool does not actually cut tubing.

Click [HERE](#) for ordering information and pictures of item # 46050-01.



Click [HERE](#) for tubing tools ordering information and pictures.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-18-2024