



Do you have a tear off crimp cap that is also magnetic – FAQ

No, we do not offer a magnetic tear off cap at this time because it requires aluminum (a soft metal) to allow easy tearing of the cap and aluminum is not magnetic. The magnetic material used in crimp caps is not soft enough to tear.

Click [HERE](#) for crimp cap ordering information and pictures.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com