

Do you have autosampler vial screw caps without septa – FAQ

Yes! We offer two types of caps without septa:

The first type is a completely closed cap where the area that would normally have a hole for the septum is covered. These are offered with **liners** instead of septa. Please see catalog numbers [9502S-09-SR](#) and [9502S-09-SR-K](#) for examples.



The second type is simply the polypropylene cap in which the septa is **Not present**. The caps are normally referred to as open hole caps. An example catalog number for this type of cap is [9502S-00-B-M](#).



Printed from the Chrom Resource Center

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Date: 05-18-2024