



Can I use vial 95025-VN-05C with crimp cap  
95025-08-1S in an oven at 105C for approximately 1  
month - FAQ

Yes, the product should performing well at this temperature and time period.

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 05-18-2024