



## What are the dead volumes of the MRQ and Fused Insert autosampler vials – FAQ

**For both of these autosampler vials, the general residual volume is about 2  $\mu$ L.**

*The dead volume is analogous to **residual volume**.*

Printed from the Chrom Resource Center

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Date: 02-05-2024