


Vials with Inserts Pre-assembled and Fused-in for Safety



Streamline your work processes and save time and money with Fused-Insert Vials instead of using separate "drop-in the vial" inserts. Fused inserts are very economical as there is no breakage or tedious filling of "hand-held" single inserts that must be carefully placed into a vial. Fused inserts are stationary, enabling you to hold the vial without damaging the insert and your sample.

Fused inserts seal directly against the septa in the cap for minimal leaking or evaporation. Wide tip inserts make it easy for rapid filling without air bubbles being formed especially when using robotics. RSA[™] glass vials and inserts will give you further confidence that your samples will not be affected by changes in pH, hydrolyze or be lost and adsorbed to the glass wall, as occurs with ordinary vials. The RSA[™] glass surface is not only the cleanest in the industry, it also has the lowest metal content ever produced for non-coated autosampler vials.

Meets or exceeds all instrument manufacturer's autosampler specifications for dimensional fit and height. Special packaging does not out-gas or add any extractables to the vials and the vials are not touched by human hands for further sample integrity.

Cat. No. Description

For more information visit www.mtc-usa.com

9532S-WCV-WT-RS Vials, 300ul clear fused insert (12 x 32 mm) screw top vials with write on patch made with RSA[™] glass, wide tip. RSA[™] Brand. 100/pack.

9532S-OCV-RS Vials, 300ul clear fused insert (12 x 32 mm) screw top vials with write on patch made with RSA[™] glass. RSA[™] Brand. 100/pack.



9158 Industrial Blvd NE Leland, NC 28451 p: 1.732.380.8900 f: 1.910.769.9435